Form - Binder

ART.NO 50209, 50220, 50204

Elastic binder which provides relief and efficient support for wounds after thoracic and abdominal surgery. Adjustable Velcro front closure makes the product easy to put on and adjust to requested compression. Can be closed right or left side. Form provides a consistent compression around the whole abdomen area and does not affect the patient's ability to breathe comfortably. The product is also suitable to be used as a stoma binder, although consider using a larger size in order to avoid too much compression. Available in black and grey in three different heights 15, 20 och 25 cm.

Level of compression: Medium, 2 out of 3.

Indications: Stabilizing girdle for stoma/parastomal hernia. Hernia prevention. Girdle/ fixation used for mobilization after surgery and as compression garment. **Contraindication:** Incarcerated hernia.

Material: 78% Polyester, 11% Polyamide, 11% Elastane. Latex free.

Washing instructions: Machine wash at 40°C (106°F). Seal the Velcros. Use washing bag. Do not bleach. Do not dry clean. Do not iron. Do not tumble dry.

Size	Measure- ment cm	Art.no Form 15 cm Black	Art.no Form 20 cm Black	Art.no Form 25 cm Black
S	70-85	50209-00-40	50220-00-40	50204-00-40
м	80-95	50209-00-50	50220-00-50	50204-00-50
L	90-105	50209-00-60	50220-00-60	50204-00-60
XL	100-115	50209-00-70	50220-00-70	50204-00-70
2XL	110-125	50209-00-80	50220-00-80	50204-00-80

Size	Measure- ment cm	Art.no Form 15 cm Grey	Art.no Form 20 cm Grey	Art.no Form 25 cm Grey
S	70-85	50209-06-40	50220-06-40	50204-06-40
М	80-95	50209-06-50	50220-06-50	50204-06-50
L	90-105	50209-06-60	50220-06-60	50204-06-60
XL	100-115	50209-06-70	50220-06-70	50204-06-70
2XL	110-125	50209-06-80	50220-06-80	50204-06-80



50204: 25 cm



50220: 20 cm



50209: 15 cm





MD Class I CE

NordiCare Ortopedi+Rehab AB Solrosvägen 1 • S-263 62 Viken, Sweden • +46 (0)42 35 22 20 • info@nordicare.se • www.nordicare.se