

Jill – Compression bolero

NEW

Seamless 4-way stretch

Flexible fit in bust
and abdomen area

Adjustable hook and
eye closure in 3 steps



NordiCare[®] +


NORDICARE ORTOPEDI + REHAB AB • Solrosvägen 1 • SE-263 62 Viken, Sverige
Tel: 042-35 22 20 • Fax: 042-35 22 21 • info@nordicare.se • www.nordicare.se

Jill – Compression bolero

NEW

Effective compression to promote fast and comfortable healing

- Recommended after arm and back surgery and liposuction.
- Compression and support for reduced wound swelling and pain relief.
- Easy application and adjustment with 3-step hook and eye closure in the front.
- Fabric in knitted seamless 4-way stretch.
- Fits different body shapes and bra sizes.
- Required uniform compression exactly where needed, targeted compression (class 1).




3-step hook & eye closure



Fully openable front



Shapes to the body



Targeted compression

Jill, black

Size	Under bust	Arm (biceps)	Waist	Art. no.
XS	63-70 cm	22-25 cm	62-68 cm	51205-00-30
S	70-77 cm	25-28 cm	68-73 cm	51205-00-40
M	75-82 cm	28-31 cm	73-79 cm	51205-00-50
L	80-87 cm	30-33 cm	79-89 cm	51205-00-60
XL	85-92 cm	31-34 cm	89-99 cm	51205-00-70
2XL	90-97 cm	33-36 cm	99-109 cm	51205-00-80



Compression class 1



Adjustable with hook & eye



Seamless



4-way stretch



Latex free

