

Seamless knitted compression

NEW

Finally in stock, place your order now!

Next generation compression garments

NordiCare's line of seamless knitted compression garments is a unique product assortment combining the most comfortable fabric together with high quality compression and functional design.



Seamless



4-way stretch



Compression class 1



40-60°C washing



NordiCare[®] 

NORDICARE ORTOPEDI + REHAB AB • Solrosvägen 1 • SE-263 62 Viken, Sweden
Phone: +46 (0)42-35 22 20 • Fax: +46 (0)42-35 22 21 • info@nordicare.se • www.nordicare.se

Seamless knitted compression

Finally in stock, place your order now!

Patsy Art.no 51201



Annie Art.no 51202



Alison Art.no 51203




Daisy Art.no 51204




Full length pants

- ✓ Hook and eye closure in both sides.
- ✓ Open crotch
- ✓ Silicone waistband.




Below knee pull-on pant

- ✓ Can be used as "Step 2 product".
- ✓ Silicone waistband.




High back body

- ✓ Hook and eye closure in both sides.
- ✓ Open crotch
- ✓ Adjustable shoulder straps.



Fully openable body

- ✓ Double material over the abdomen/pocket.
- ✓ Hook and eye closure in both sides and crotch.
- ✓ Silicone waistband.



Patsy, Annie, Alison, Daisy

	Size	Waist cm	Hip cm	Thigh cm
	X Small	62-68	86-91	46-50
	Small	69-75	92-97	50-54
	Medium	76-83	98-103	54-58
	Large	84-92	104-111	58-62
	X Large	93-101	112-119	62-66
	2X Large	102-110	120-127	66-70